

Procedure: Knee Arthroscopy & partial menisectomy +/- Meniscal Repair	OPCS 4.3 Code: W8200 +/- W8300
--	---

<p style="text-align: center;">Sets:</p> <p>Arthroscopy Set. 30 degree 5-5.5 mm dia Smith & Nephew scope.</p> <p>Have available but not open:</p> <ul style="list-style-type: none"> • Arthroscopy Punches • Smith & Nephew power arthroscopic shaver and foot control 	<p style="text-align: center;">Extras:</p> <p>Have available but not open:</p> <ul style="list-style-type: none"> • Incisor shaver blades (Smith & Nephew) <ul style="list-style-type: none"> - 3.5 incisor -7205312 - 4.5 incisor - 7205345 • Arthrowand - <ul style="list-style-type: none"> - Super Multivac 50 (Arthrocare) - AS4630-01 • Disposafe Plus Sharps mat. (Amba Medical Ref. QF-1001) • 20ml syringe • 21G Green needle • Suction Tubing • Y-giving Set • Smith & Nephew 5x Curved FAST-FIX Meniscal Repair System REF 7207877 • Smith & Nephew 1x FAST-FIX Meniscal Knot Pusher and Suture Cutter REF 7209084 	<p style="text-align: center;">Sutures:</p> <ul style="list-style-type: none"> • 3/0 Nylon – W320
<p style="text-align: center;">Swabs:</p> <ul style="list-style-type: none"> • Small x 1 		<p style="text-align: center;">Pharmacy:</p> <ul style="list-style-type: none"> • 20ml 0.5% Bupivacaine with 1:200,000 Adrenaline • Warm 3L Saline
<p style="text-align: center;">Blades:</p> <ul style="list-style-type: none"> • 11 x 1 • Blade handle 		

Prep: Alcoholic Betadine

<p style="text-align: center;">Drape:</p> <ul style="list-style-type: none"> • U Shaped Drape with adhesive aperture PE Drape P3 Medical Ref: FXS01ORT0001 • Camera Drape (Amba, Clinicon Ref. ECN 150-250) • Arthroscopy drape w/pouch Molnycke • Adhesive edge Utility Drape for foot • 6" Crepe bandage to wrap utility 	<p style="text-align: center;">Dressings:</p> <ul style="list-style-type: none"> • 2x Small Primapore 9 x 15 • 6" Velband • 6" Crepe Bandage • Elastoplast
<p style="text-align: center;">Table Set-up:</p> <ul style="list-style-type: none"> • Supine with lateral support level with lower edge of tourniquet and by patient's thigh (Amatech Supports) • Heel gel support contralateral limb 	<p style="text-align: center;">Theatre Set-up:</p> <ul style="list-style-type: none"> • Razor for shaving leg & Yellow bin bag • Tourniquet 350 mmHG– Velband underneath • Rees-Davis exsanguinator • Intellijet Stack with Arthrocare machine • Smith & Nephew Camera stack system

General Information:

- Always recorded for the patient in video and stills
- Input Mr. J. Hardy M.D for still image labels
- Clear video and images files for next patient after recording
- Medical Students and guests observe from outside of laminar flow enclosure